


SQUIDTM peri



The New Liquid Embolic Device

SQUIDPERI is a non adhesive liquid embolic agent for embolization of Peripheral Arteriovenous Malformations (AVM). It is composed of EVOH (Ethylene Vinyl Alcohol Copolymer) with suspended micronized Tantalum powder for radiopacity, and DMSO (Dimethyl Sulfoxide) solvent. SQUIDPERI must be injected through a compatible microcatheter.

SQUIDPERI EMBOLIZATION RANGE:
the POWER OF VISIBILITY, the CHOICE OF FLUIDNESS



Route des Avouillons 30 | 1196 Gland | SWITZERLAND
Tel: +41 26 321 44 80 | info@emboflu.ch | www.emboflu.ch

A UNIQUE MICRONIZATION PROCESS

EMBOFLU developed a specific micronization process to minimize Tantalum powder grain size in the SQUIDPERI suspension. This enables a SLOWER precipitation of the radiopaque powder which stays LONGER in the SQUIDPERI suspension.

**Micronized grain size
of Tantalum powder
(SQUIDPERI 18)**

Gauss-like repartition
of the Tantalum powder
still in suspension

SQUIDPERI →
Tantalum powder
deposit



**Standard grain size
of Tantalum powder**

← Double the thickness
of SQUIDPERI

Image taken 15 mins after shaking

BENEFITS of the SQUIDPERI unique micronization process

- **High homogeneity of SQUIDPERI suspension** - eliminates aggregates formation which can cause microcatheter blockage.
- **HOMOGENIC radiopacity** - reduces discrepancy between saturated radiopaque zones and non-saturated zones, for improved assessment of AVM angioarchitecture.
- **High STABILITY over time** - allows improved VISIBILITY for longer INJECTION TIMES


4 SQUIDPERI FORMULAS

STANDARD VISCOSITY

SQUIDPERI 18: Standard version of SQUIDPERI for standard AVM embolizations.


SQUIDPERI 18LD (Low Density): 30% less Tantalum than standard SQUIDPERI 18.

- For a better assessment of the AVM vasculature and the volume of liquid embolic injected.
- To avoid the over saturated radiopaque injected zones (flash-effect).



AVM of the hip.

Courtesy of Dr. Patrick Brouwer, Rotterdam ,Holland



pre-surgery injection of SQUID 12 into a tumor.

Courtesy of Dr. Patrick Brouwer, Rotterdam ,Holland

LOW VISCOSITY

SQUIDPERI 12: With lower viscosity, this version is more fluid than the standard formula and allows:

- Deeper penetration into the nidus.
- Reaching distal microvessels and injection through small feeders.

SQUIDPERI 12LD (Low Density):

This formula has the same viscosity as SQUIDPERI 12, but 30% less Tantalum.

- To enable assessment of the AVM vasculature and the volume of liquid embolic injected.



SQUIDTM peri

The New Liquid Embolic Device

Ordering Information

References	Description	Content of each reference
SQUIDPERI 18	standard viscosity SQUIDPERI	● one 1.5ml vial of SQUIDPERI
SQUIDPERI 18LD	standard viscosity and low density SQUIDPERI: 30% less Tantalum than SQUIDPERI 18	● one 1.5 ml vial of DMSO ● one 1cc Blue syringe for DMSO
SQUIDPERI 12	low viscosity SQUIDPERI	● two 1cc White syringes for SQUIDPERI
SQUIDPERI 12LD	low viscosity and low density SQUIDPERI: 30% less Tantalum than SQUIDPERI 12	● two syringe adapters



SQUIDPERI EMBOLIZATION RANGE:
the POWER OF VISIBILITY, the CHOICE OF FLUIDNESS

