

introducer sheath (long)

ballast

Reference	Usable length (cm)	Proximal O.D.	Distal O.D.	I.D.	Hydrophilic coating (cm)
BALLAST80	80	0.106"	0.100"	0.088"	20
BALLAST90	90				
BALLAST100	100				

IVA

Reference	Usable length (cm)	I.D.
IVA3F	7	3F
IVA3F50	50	3F
IVA5F45	45	5F
IVA6F80	80	6F
IVA6F80ST	80	6F
IVA6F80ST-MP	Supple Tip = 50 mm	6F



guiding catheter

fargomini

Reference	Proximal O.D. - Distal O.D.	I.D.	Total length (cm)	Super-supple distal tip length (cm)	Curve
FRGMIN4,2F120	4,2 F to 3,9 F	1,02 mm (.040")	120	2.5	Straight (D)/ Multipurpose (MP)
FRGMIN4,2F120MP			135	15	
FRGMIN4,2F135					
FRGMIN4,2F135MP					

fargo

Reference	Proximal O.D. - Distal O.D.	I.D.	Total length (cm)	Super-supple distal tip length (cm)	Curve
FRG6F105_8	6,0 F to 4,9 F	1,46 mm (.058")	105	8	Straight (D)/ Multipurpose (MP)
FRG6F105_8MP			115		
FRG6F115_8					
FRG6F115_8MP					
FRG6F125_8			125		
FRG6F125_8MP			135	15	
FRG6F135_15					
FRG6F135_15MP					

fargomax

Reference	Proximal O.D. - Distal O.D.	I.D.	Total length (cm)	Super-supple distal tip length (cm)	Curve
FRGMAX6F95_8	6,0F to 6,0F	1,78 mm (.070")	95	8	Straight (D)/ Multipurpose (MP)
FRGMAX6F95_8MP			105		
FRGMAX6F105_8					
FRGMAX6F105_8MP					
FRGMAX6F115_8			115		
FRGMAX6F115_8MP			125		
FRGMAX6F125_8					
FRGMAX6F125_8MP					

