


progressive suppleness microcatheter

magic

Reference	Total length (cm)	Guidewire max	Distal tip O.D.	Turquoise pursil O.D. (F) x L (cm)	White pursil O.D. (F) x L (cm)	Pink pursil O.D. (F) x L (cm)	Torver O.D. (F) x L (cm)
MAGIC1,2F	165	.008"	1,2F	1,2 x 12	1,8 x 3	2,4 x 30	2,7 x 114
MAGIC1,2FM	165		1,2F	1,2 x 3	1,8 x 12	2,4 x 30	
MAGIC1,2FM180	180		1,2F	1,2 x 3	1,8 x 17	2,4 x 40	
MA1,2FOLIVE	165		1,2F	1,2 x 12	1,8 x 3	2,4 x 30	
MAGIC1,5F	160	.009"	1,5F	/	1,5 x 15	2,4 x 25	2,7 x 114
MAGIC1,5MP	165		1,5F	/	1,5 x 20		
MAGICSTD	155	.010"	1,8F	/	1,8 x 10	2,4 x 25	2,7 x 114
MAGICMP	165		1,8F	/	1,8 x 20		

baltacci

Reference	Total length (cm)	Guidewire max	Super supple part O.D. (F) x L (cm)	Supple part length (cm)	Torver O.D. (F) x L (cm)
BALTACCI1,2F	165	.007"	1,2 x 8	22	2,7 x 129
BALTACCI1,2FP	100				2,7 x 64
BALTACCI1,5F	165	.009"	1,5 x 10	20	2,7 x 129
BALTACCI1,8F	165	.009"	1,8 x 10	20	2,7 x 129
BALTACCI1,8FP	100			14	2,7 x 70



detachable tip microcatheter

sonic

Reference	Total length (cm)	Distal tip O.D.	Distal tip I.D. (mm)	Fusecath length (cm)	Dead space		Max reflux Squid/ Onyx® (cm)	Included guidewire
					with adapter (ml)	without adapter (ml)		
SONIC1,2F15	167	1,2F	0.23	1.5	0.24	0.28	2	.007"
SONIC1,2F25	167	1,2F	0.23	2.5	0.25	0.29	3	.007"
SONIC1,2F35	167	1,2F	0.23	3.5	0.25	0.29	4	.007"
SONIC1,2F25.190	190	1,2F	0.23	2.5	0.28	0.32	3	.007"
SONIC1,5F15	167	1,5F	0.3	1.5	0.26	0.3	2	.008"
SONIC1,5F25	167	1,5F	0.3	2.5	0.26	0.3	3	.008"
SONIC1,5F35	167	1,5F	0.3	3.5	0.25	0.29	4	.008"
SONIC1,5F45	167	1,5F	0.3	4.5	0.26	0.3	5	.008"

