


venous remodelling balloon

copernic RC

Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Catheter outer maximum diameter	Guidewire recommended	Compatible guiding catheter
COPERNIC8X80RC	160	80	8	4,35F (1,45mm)	compatible with .014" guidewire	Fargomax 6F (ID: 1,78 mm / .070")
COPERNIC10X80RC			10	5F (1,6mm)		



detachable balloon

goldballoon

Reference	Volume (mL)	Inflated DxL (mm)	Collar Ø (mm)	Minimum guiding catheter internal Ø (mm)
GOLDBAL1	0.25	7x10	0.9	1.75
GOLDBAL2	0.6	7x20		1.75
GOLDBAL3	0.5	9x11		2.2
GOLDBAL4	0.75	9x16		2.2
GOLDBAL5	2.5	12x28		2.8
GOLDBAL6	3	15x25	1.05	3

catheters for detachable balloons

magic

Reference	ø Proximal / distal (F)	Length (cm)
MABDTE	2,7 / 1,8	165

baltacci

Reference	ø Proximal / distal (F)	Length (cm)	Max ø guidewire
BALTACCIBDPE	2,7 / 2,4	163	.009"
BALTACCIBDPE100	2,7 / 2,5	94	.009"

