

intracranial self-expandable flow diverter

silk vista baby


| Reference | Vessel Ø (mm) | Stent's length at nominal Ø (mm) | Unconstrained stent | | Delivery catheters |
|----------------|---------------|----------------------------------|---------------------|-------------|--------------------|
| | | | Ø (mm) | length (mm) | |
| SILK_V_2,25x10 | 2.25 | 10.5 | 2.5 | 8 | .017" |
| SILK_V_2,25x15 | | 16 | | 12 | |
| SILK_V_2,25x20 | | 22 | | 15.5 | |
| SILK_V_2,75x10 | 2.75 | 12.5 | 3 | 9 | .017" |
| SILK_V_2,75x15 | | 17 | | 12 | |
| SILK_V_2,75x20 | | 22 | | 15 | |
| SILK_V_2,75x25 | | 26.5 | | 18.5 | |
| SILK_V_3,25x10 | 3.25 | 11 | 3.5 | 8.5 | .017" |
| SILK_V_3,25x15 | | 16.5 | | 12.5 | |
| SILK_V_3,25x20 | | 21 | | 15 | |
| SILK_V_3,25x25 | | 26 | | 18 | |

silk+

| Reference | Vessel Ø (mm) | Stent's length at nominal Ø (mm) | Unconstrained stent | | Delivery catheters |
|-------------|---------------|----------------------------------|---------------------|-------------|--------------------|
| | | | Ø (mm) | length (mm) | |
| SILK2,0x10* | 1,50 to 2,25 | 10 | 2.5 | 6 | vasco+21 |
| SILK2,0x15 | | 15 | | 10 | |
| SILK2,0x20 | | 20 | | 11 | |
| SILK2,5x10* | 2,00 to 2,75 | 10 | 3 | 7 | vasco+21 |
| SILK2,5x15 | | 15 | | 10 | |
| SILK2,5x20 | | 20 | | 13 | |
| SILK2,5x25 | | 25 | | 16 | |
| SILK3,0x10* | 2,50 to 3,25 | 10 | 3.5 | 7 | vasco+21 |
| SILK3,0x15 | | 15 | | 10 | |
| SILK3,0x20 | | 20 | | 13 | |
| SILK3,0x25 | | 25 | | 16 | |
| SILK3,0x30 | | 30 | | 19 | |
| SILK3,5x10* | 3,00 to 3,75 | 10 | 4 | 9 | vasco+21 |
| SILK3,5x15 | | 15 | | 11 | |
| SILK3,5x20 | | 20 | | 12 | |
| SILK3,5x25 | | 25 | | 14 | |
| SILK3,5x30 | | 30 | | 15 | |
| SILK3,5x35 | | 35 | | 17 | |
| SILK4,0x15 | 3,50 to 4,25 | 15 | 4.5 | 10 | vasco+21 |
| SILK4,0x20 | | 20 | | 12 | |
| SILK4,0x25 | | 25 | | 14 | |
| SILK4,0x30 | | 30 | | 16 | |
| SILK4,0x35 | | 35 | | 18 | |
| SILK4,0x40 | | 40 | 20 | | |
| SILK4,5x15 | 4,00 to 4,75 | 15 | 5 | 8 | vasco+21 |
| SILK4,5x20 | | 20 | | 10 | |
| SILK4,5x25 | | 25 | | 12 | |
| SILK4,5x30 | | 30 | | 15 | |
| SILK4,5x35 | | 35 | | 17 | |
| SILK4,5x40 | | 40 | | 19 | |
| SILK5,0x20* | 4,50 to 5,25 | 20 | 5.5 | 9 | vasco+25 |
| SILK5,0x25 | | 25 | | 11 | |
| SILK5,0x30 | | 30 | | 13 | |
| SILK5,0x40 | | 40 | | 17 | |
| SILK5,5x25 | 5,00 to 5,75 | 25 | 6 | 11 | vasco+25 |
| SILK5,5x30 | | 30 | | 13 | |
| SILK5,5x40 | | 40 | | 17 | |
| SILK5,5x50* | | 50 | | 23 | |

tapered silk+

| Reference | Vessel Ø (mm) | Vessel Ø distal (mm) | Stent's length at nominal Ø proximal (mm) | Stent length at nominal Ø distal (mm) | Unconstrained stent | |
|-----------------|---------------|----------------------|---|---------------------------------------|---------------------|-------------|
| | | | | | Ø (mm) | length (mm) |
| SILKP4,0D3,0x30 | 3,50 to 4,25 | 2,50 to 3,25 | 12 | 10 | 4,5/3,5 | 15 |
| SILKP4,5D3,0x25 | 4,00 to 4,75 | 2,50 to 3,25 | 12 | 10 | 5,0/3,5 | 20 |
| SILKP4,5D3,5x30 | 4,00 to 4,75 | 3,00 to 3,75 | 11 | 15 | 5,0/4,0 | 16 |



*manufactured on demand